

Processing Breath-by-Breath Data From Indirect Calorimetry

Robert A. Robergs, Ph.D., FASEP, EPC

Director: Exercise Physiology Laboratories, Exercise Science Program, Department of Health, Exercise and Sports Science, The University of New Mexico, Albuquerque, NM, USA. Phone: (505) 277-2658; FAX: (505) 277-9742; Email: rrobergs@unm.edu

Nutrition, Exercise and Sports Information Services (NEXSIS); www.NEXSIS.org

Overview

The program you purchased and are about to use was written with LabVIEW, a data flow programming language from National Instruments (Austin, Texas). A LabVIEW program is referred to as a virtual instrument (VI). This program has one main control program, and within this are several other sub-programs, or subVIs.

With this program you will be able to process breath-by-breath data using either of user-selected time or breath averages, as well as a low frequency digital filter at user-selected low frequency (Hz) cut-off values. Figure 0 presents a photograph of a subject performing an incremental exercise test using breath-by-breath technology. The program also supports the selection of regions within the data set to apply segmental linear regression to more objectively identify threshold changes, such as for the ventilation or lactate threshold. Data sets from within the imported text file can be saved after processing a text file.

LabVIEW programming is based on what are called the Front Panel and the Block Diagram. The Block Diagram is for the programmer for developing the program. The Front Panel is essentially the user interface.



Figure 0. A female endurance trained subject completing a test of VO_2 max in the Exercise Physiology Laboratories of the University of New Mexico.

Figure 1 provides an image of the Front Panel for this program.

Note that there are 8 options, and since you are reading this, you have managed to successfully migrate to this Program Help window. Alternatively, you may have gone straight to the pdf file that is called by the “Program Help” button after you installed the software. Note that the options are grouped into String File, Data File, Header File, Process Data and Thresholds sections, and all are color coded. The Program Help icon is separate. When you click on an option, you will start a subVI with a specific and different Front Panel to this “Home” window. These features are explained below.

Program Instructions

As you now know, the Program Instructions selection takes you to a window that displays an Acrobat Reader display containing a pdf file of this instruction document. Scroll down and up through the pdf document with this Acrobat Reader display and learn how to use and get the most from this program. If you do not see this pdf when you run this program, you do not have the correct version of Adobe Acrobat Reader. My experience with this program and different versions of Adobe Acrobat Reader has revealed that this program works best with Adobe Acrobat Reader version 7.

String Files

A string file is a text file that contains words, or both words and numbers. For LabVIEW to read a file with words in it, additional programming needs to occur in order to access numbers vs. word text. Hence the distinction between string and data files within this program. Most often the data file you import that contains words will have header names, with additional information. As you will see, you need to edit this raw file so that it contains only 1 row of header words (1 header per column) and the remaining content being all numbers.




Figure 1. An image of the Front Panel window of the Program and the initial program selection features.

Get String File

Figure 2 reveals the Front Diagram of the “Get String File” subVI, along with the file selection window that appears when the sub-VI is first called.

After a string file is selected, the text and data appear in the array as shown in Figure 3. This data is from a typical text file generated by a commercial system for indirect calorimetry. It is recommended that you check all the columns and rows of this data file by scrolling through the row and column scroll bars to the top left of the data array. This will reveal what columns and rows that you need to delete in the next subVI.

In the data file of Figure 3, the initial text and data of columns 0 to 8 need deleting, as do the contents of row 1 and 2. Remember number counting starts at 0 in LabVIEW and all programming. When this is done, as shown in Figure 4, the Time data column is the 0 column, with the header “t” at row 0, column 0.

Click “Continue” to go back to the control VI so you can select the “Edit String File” subVI and complete the necessary editing.




Figure 2. An image of the Front Panel of the “Get String File” subVI.




Figure 3. Typical string data imported into the program.

Edit String File

Figure 4 reveals the Front Diagram of the “Edit String File” subVI. The text and data of Figure 4 reveal the single header per column of selected, as explained above. Note that the scroll bars for rows and columns are still on “0”. As there are no additional rows that need deleting at the top of the file, clicking on the row scroll bar would increment rows one at a time until you get to the bottom/end of the data file. Alternatively, highlighting the row number of the scroll bar and typing a given row number and clicking off the scroll bar would immediately increment the top row to this row number. The same is true for the column scroll bar.

Delete Data

The “Delete Data” window can then be used to select a row or column to delete. Clicking either of the green circular icons causes the deletion. It is important to note that the selected row or column number does not change after each deletion. This way, you can select the first row or column of a block of data that needs deleting, and simply continue to click on the green circular icon until all rows or columns have been deleted.

Undo Last Delete

The “Undo Last Delete” can only do what it says, undo the LAST deletion. You can then insert this column or row where you want by selecting a given row or column. This feature also lets you rearrange the data set, by row or column. The bottom line though is to not delete more than one row or column at a time before you realize that you have made a mistake. If you do this, you will have to start the program again to re-edit your data file, as you will not be able to retrieve rows or columns deleted prior to the last delete. Sorry!!!

Edit String File Click to Save & Continue Continue

Data	t	Rf	VT	VE	VO2	VCO2	O2exp	CO2exp
0	0.000000	19.60784314	2.81308392	55.15850824	2436.002895	2315.75031	427.3193855	154.2757667
0	0.050000	26.78571429	1.588673074	42.55374304	1877.576866	1747.188046	238.8621494	87.95351249
0	0.066667	35.71428571	1.546839036	55.2442513	2747.995647	2327.796507	223.4549795	87.74944402
0	0.100000	37.5	1.425418294	53.45318603	2573.708299	2146.52303	207.741707	78.15822572
0	0.150000	22.90076336	2.536571138	58.08941536	2693.895408	2457.008864	376.0981984	143.1540348
0	0.183333	29.12621359	1.890694416	55.06876939	2553.376075	2345.964397	280.7982209	107.4420518
0	0.216667	21.50537634	1.521330477	32.71678445	1844.410033	1510.00858	207.1295016	94.68777215
0	0.266667	26.43171806	3.086535676	81.58244078	4102.641818	3704.05287	445.0733428	184.2738324
0	0.300000	30.15075377	2.528408399	76.23341905	3364.239123	3268.473733	382.6283896	142.0316582
0	0.316667	34.88372093	1.742744772	60.79342227	2757.02719	2633.11045	259.2689968	101.0138949
0	0.350000	35.92814371	1.944772562	69.87206808	3022.181416	2942.952931	293.348432	110.9112159
0	0.383333	39.47368421	1.964159067	77.53259474	3411.762947	3346.902743	294.8789456	114.3803799
0	0.400000	34.09090909	2.084559467	71.06452727	2961.008647	3028.520026	317.3264778	120.8085369
0	0.433333	36.58536585	2.174349595	79.54937544	3361.546672	3386.44526	331.611271	124.0736325
0	0.466667	38.70967742	2.161085145	83.65490882	3252.920414	3583.771757	334.3661954	126.8285569

Note: Deleting rows or columns, and then using the "insert where" Undo icon also allows you to re-arrange the spreadsheet rows and columns.

Delete Data

What row do you want to delete? ● Row ● click to delete

(remember LabVIEW numbering starts at 0)

What column do you want to delete? ● Column ● click to delete

Make sure there is 1 row header for each data column

Undo Last Delete

● Undo

insert where?

● Undo

Convert Time to Min Fraction

Min.sec

For time stamp data only
(00:00:00)
(-H:MM:SS)

Figure 4. Data of Figure 3 that has been edited to show column headers in row 0 and data in all remaining rows.

Convert Time to Minute Fraction

The “Convert Time to Minute Fraction” option is for data files that have time stamp data for the time column. Look back to the time (t) data of Figure 3 (last column on right). This is time stamp data. This option will convert time stamp to minute fraction data, as you see in the first column (column 0) of Figure 3. **For time data that is not already in a minute fraction, or a second fraction, or time stamp, you will need to convert this time data column to a time (minute or second) fraction in Excel prior to using the text file in this program.**

Clicking the “Continue” icon first launches a prompt to name and save the edited string file a destination path. Once this occurs, the program returns to the control front panel.

Data Files

Some VO₂ systems save only data to a text file. Here you simply import this data file and then after editing it, you make your own headers for each column.

Get Data File

Figure 5 reveals the Front Diagram of the “Get Data File” subVI. This subVI is similar to the “Get String File” subVI, as revealed in Figures 2 and 3, only it is for data text files only. In this example, only 5 columns are shown (0-4) in the data set, and these are time, expired ventilation, VO₂ (L/min), VO₂ (mL/kg/min) and heart rate.

Edit Data File

Figure 6 reveals the Front Panel of the “Edit Data File” subVI. This is similar to the “Edit String File” of Figure 4.

The only difference with this “Edit Data File” subVI is that the Continue button also saves a data text file of the edited data just in case you need this edited file for other reasons, such as for importing into a graphics program. Once you click on this button icon you will be prompted to name this text file and select the path where you want




Figure 5. The data retrieved from a text file containing 5 columns of data.

to save it to. As indicated, you should only be working with numbers, and as such, you can only save this data to a file if it contains only numbers.

Name Columns

The “Name Columns” subVI is for work with data files only. As string files have their own text headers, and these are saved in the programming when editing the string file, ***you do not need to use this subVI if you worked with a string file.***

When using data files, you need to create headers, or column names, to be able to use the data processing subVI. You do this with this subVI, as shown in Figure 7. Note that there are “Enter” and “Edit” selections. To name a column, select a column number and type the name in the Column Name entry box. When you click the “Enter” icon, you will see the column number and name in the Header 2D array. In Figure 7 you will see that I named the first three columns (0-2), and then edited the time header to show the units of time. As the column number does not increment after each entry, make sure you remember to change the column with each new column header.


To save the column names (headers) to a file for latter use in this program when working with a different data set from the same system and final data column features, click on the “Save Header File” button. This will prompt you to enter a name and select the path where you want this file saved. If you edit your VO₂ data the same, then I would use this saved header file rather than use the “Name Columns” subVI each time I processed data. This will save you time when processing multiple data files.

The screenshot shows the 'Edit Data File' subVI interface. At the top, it displays 'Array 2' with a 2D array of numerical data. Below the array, there are controls for deleting data. The 'Delete Data' section has two input fields: 'What row do you want to delete?' and 'What column do you want to delete?'. Both fields have a value of 0. There are also 'Row' and 'Column' labels and a 'click to delete' button. A note below these fields says 'Make sure only numbers remain'. The 'Undo Last Delete' section has two input fields: 'insert deleted column where?' and another field with a value of 0. There are also 'Edit' buttons for both sections. A 'Continue' button is located on the right side of the interface.

0.1	9.474	14.003	9.536	14.094	0.273	0.212	3.389	2.628	0.775	34.683	44.724	68.458	0.182
0.156	16.347	24.161	16.429	24.282	0.437	0.355	5.425	4.408	0.813	37.389	46.015	64.388	0.184
0.204	15.346	22.682	15.412	22.778	0.466	0.401	5.785	4.975	0.86	32.912	38.27	65.368	0.18
0.263	16.176	23.908	16.408	24.251	0.756	0.524	9.382	6.504	0.693	21.393	30.858	75.1	0.166
0.318	10.054	14.86	10.187	15.057	0.379	0.246	4.703	3.053	0.649	26.526	40.855	76.517	0.174
0.376	10.856	16.045	10.881	16.083	0.268	0.243	3.325	3.012	0.906	40.514	44.721	67.458	0.185
0.431	10.469	15.474	10.535	15.571	0.28	0.214	3.472	2.658	0.765	37.415	48.877	64.413	0.184
0.491	13.17	19.465	13.202	19.513	0.363	0.33	4.501	4.099	0.911	36.3	39.867	67.387	0.182
0.543	13.702	20.252	13.811	20.413	0.464	0.356	5.763	4.413	0.766	29.5	38.527	71.417	0.177
0.592	11.146	16.474	11.15	16.479	0.347	0.344	4.311	4.269	0.99	32.078	32.393	71.457	0.178
0.649	10.006	14.789	10.037	14.834	0.271	0.24	3.363	2.984	0.887	36.92	41.609	69.432	0.183
0.711	21.04	31.097	21.153	31.264	0.557	0.444	6.909	5.508	0.797	37.781	47.391	70.432	0.184
0.802	31.234	46.164	31.725	46.889	1.803	1.312	22.367	16.282	0.728	17.325	23.801	76.542	0.155
0.862	10.105	14.935	10.178	15.043	0.33	0.258	4.098	3.197	0.78	30.596	39.213	71.874	0.178
0.931	13.112	19.38	13.174	19.471	0.453	0.392	5.626	4.858	0.863	28.915	33.49	73.462	0.176
0.983	8.138	12.028	8.129	12.014	0.212	0.221	2.626	2.737	1.042	38.454	36.892	74.472	0.183

Note: Deleting rows or columns, and then undeleting also allows you to re-arrange the spreadsheet rows and columns.


Figure 6. The data retrieved from a large text file prior to any editing. As explained in the text, editing could consist of deleting columns, or rearranging the columns.



Get Header File

Figure 8 reveals the Front Panel of the "Get Header File" subVI.

When working with data files, and when using a header file from a previous data processing, this allows you to import this file and converts it to a suitable format to use in the "Process Data" subVI.



Process Data

Figure 9 reveals the Front Panel for the “Process Data” subVI, which also functions to assist in detecting maximal values.

This data is from an incremental exercise test to VO_2 max. Note the four user interactive controls in the top left for Select Processing, Select X Variable, Select Y Variable, and Breath Number. Actually, the Breath Number control changes depending on what processing is selected. In this example, a Breath average is selected, and so a control for Breath Number appears. If the selection was for a time average, a control for the time of the average appears. If a Digital filter is selected for processing, then a control for the lower Hz cut-off appears. As previously explained, the headers for the data columns are used in the Select X or Y Variable controls. If “Nothing” is selected for the processing, no data appears on the graphs. A breath average of 1 breath would be the raw data when using breath-by-breath data.




Figure 9. An image of the Front Panel of the “Process Data” subVI.

The “Max” number above the right graph reveals the highest data point in the data set after processing. My argument is that the highest data point after your selected processing would be the peak data point, or either VO_2 peak or VO_2 max when processing VO_2 data. However, you can use whatever method you want for detecting peak values. Just be aware that when using time averages, even a 30 s time average severely over-processes the data and falsely lowers maximal values. This will be self evident when you play with the program. You will see that a 0.08 Hz digital filter, a 7 breath average and a 15 s average give similar processing.

The “Save to Data File” icon can be used to save each time and Y variable data set to a text file after it is processed.

Thresholds

Figure 10 presents the front panel image of the “Thresholds” option. Note the top figure for first processing the data, which is identical to the options of Figure 9. However, in this graph, there are also 6 cursor bars, with a cursor legends bar underneath the graph on the right side. The 6 cursors actually represent 3 pairs of cursors, as revealed by the cursor names, each one having a start and end reference (e.g. start0 and end0). To remove the cursors from the graph, simply enter data values for X and Y coordinates that are not on the visible region of the X and Y axes of the graph. For example, as all cursors have Y values set to 200, which is off the Y scale of the graph, no horizontal lines are visible on the graph.

In the data example of Figure 10 I have used the cursors to assess any slope change in VCO_2 . I prefer to use $VEVO_2$ and $VEVCO_2$ to detect the ventilation threshold. I start with the 3rd cursor pair and position these at the end of the data set for the last segment of data. I then position the 2nd cursor pair over the prior segment, and the initial cursor pair across what is the next segment. The order of the cursors is important, and always order from left to right.

The graph on the bottom left is the raw VO_2 data, and it used to generate a linear regression between the region bordered by the blue and yellow cursor. Note that you select the column that contains the VO_2 data, and remember that the first column is column 0. When the “Compute” circular icon is clicked, the segments generated on the top main graph by the cursors are used to detect two thresholds where the three lines intersect. These are “Threshold 1” and “Threshold 2”, respectively. The VO_2 at these threshold times are computed from the VO_2 -time linear regression and




Figure 10. An image of the Front Panel of the “Thresholds” subVI for expired ventilation data from a test to VO_{2max} .

presented in their respective indicators. Diff1 and Diff2 represent the differences in processing math used to detect the thresholds. They represent minimal differences between VO_2 predictions based on time to show where two linear segments are close enough to being equal, and hence intersect. I have found that the default of 0.075 s works well for most VO_2 breath-by-breath data sets.

Be careful when using this application for VO_2 data sets where the VO_2 slope of the test is not consistent across the times of the thresholds. When this happens, select a segment of the VO_2 curve that only spans single thresholds, which will require you to repeat the processing of the VO_2 region for each threshold time.

The red plateau box changes to green when the VO_2 data changes by less than 50 ml/min over the last 30 s of the data set. You need to make sure the data set ends at the end of the test (contains no recovery data) for this to be valid.